
# COD Setups COD (ISO 15705:2002) COD Photometer

#### Determination of the chemical oxygen demand index (ST-COD)

Small-scale sealed-tube Total range 0 - 15000 mg/l



The chemical oxygen demand, ST-COD value (ST = small scale **s**ealed **t**ube), of water as determined by this dichromate method can be considered as an estimate of the theoretical oxygen demand, i.e. the amount of oxygen consumed in total chemical oxidation of the organic constituents present in the water.

#### **COD Photometers**

With a measuring range from 0 to 15,000 mg/l  $O_2$ , the Lovibond® COD photometers are suitable for waste water testing.

Two LEDs light sources with long-term stability ( $\lambda_1 = 610 \text{ nm}$ ;  $\lambda_2 = 430 \text{ nm}$ , according to ISO 15705:2002), a waterproof sample chamber, a large digital display, and the user-friendly keypad ensure maximum operating reliability and convenience.

MD 100 COD (in case)	Order code: 27 61 20
MD 110 COD (in case)	Order code: 296 12 02
MD 200 COD (in case)	Order code: 289 25 02
MD 600 COD (in case)	Order code: 21 40 20
MD 610 COD (in case)	Order code: 21 40 25

#### Ranges

66

 $0 - 150 \text{ mg/l } O_2 \pm 3.5 \text{ %}^*) \text{ FS}$   $0 - 1500 \text{ mg/l } O_2 \pm 3.5 \text{ %}^*) \text{ FS}$  $0 - 15000 \text{ mg/l } O_2 \pm 3.5 \text{ %}^*) \text{ FS}$ 

\* tolerance based on the use of potassium-hydrogenephthalate standards (DIN 38409)



#### **Setups COD**

The Lovibond® COD Setups allow highly sensitive and precise water testing with minimum effort. They measure the ST-COD concentration by photometric detection employing a linear relationship between absorbance and concentration.

After adding the sample to a Lovibond® COD tube test (LR, MR according to ISO 15705:2002), it is heated in the reactor for two hours at 150 °C and then analysed in the photometer.

The COD Setups comprise the photometer, 25 tube tests for each of the two lower measuring ranges, a reactor for sample digestion and a vial stand.

#### **COD Setups**

**MD 100 COD** Order code: 27 61 30 Instrument in carrying case

**MD 110 COD** Order code: 29 61 302 Instrument in carrying case

**MD 200 COD** Order code: 289 26 02 Instrument in carrying case

**MD 600 COD** Order code: 21 40 40 Instrument in carrying case

**MD 610 COD** Order code: 21 40 41 Instrument in carrying case

#### **Delivery Content**

- adapter for round vials ø 16 mm
- 2 sets of tube tests 0-150 mg/l 0-1500 mg/l
- thermoreactor RD 125
- tube stand
- 2 syringes 1 ml, 2 ml
- batteries
- warranty information
- certificate (COC)
- instruction manual

#### **COD VARIO tube tests**

The Lovibond® COD VARIO tube tests are available for the measuring ranges 0-150 mg/l  $O_2$ , 0-1500 mg/l  $O_2$  and 0-15000 mg/l  $O_2$ .

Their chemical properties and a 16 mm tube diameter make them compatible to Hach® devices.\*

Tube tests	Order code
0-150 mg/l O₂	
(25 pc.), mercury free**	2 42 07 10
(25 pc.)	2 42 07 20
(150 pc.)	2 42 07 25
0-1500 mg/l O <sub>2</sub>	
(25 pc.), mercury free**	2 42 07 11
(150 pc.), mercury free**	2 42 07 16
(25 pc.),	2 42 07 21
(150 pc.)	2 42 07 26
0-15000 mg/l O <sub>2</sub>	
(25 pc.), mercury free**	2 42 07 12
(25 pc.)	2 42 07 22
(150 pc.)	2 42 07 27
** without chlorido romoval	

<sup>\*\*</sup> without chloride removal

#### Standard solutions

Standard solutions are solutions with a defined concentration and are provided to check the operation methods and devices of the cuvette tests as well as the condition of optical filters and the instrument.

Standard solution	Quantity	Code
100 mg/l COD	30 ml	2 42 08 03
500 mg/l COD	30 ml	2 42 08 04
5000 mg/l COD	10 ml	2 42 08 05

#### Highlights

- ST-COD sealed tubes ready for use
- Suppression of chloride interference up to 1000 mg/l (LR & MR) up to 10000 mg/l (HR)
- Mercury free tube tests, in absence of chloride interference
- 3 ranges:

Low range:

0 - 150 mg/l, meets ISO 15705:2002 Middle range:

0 - 1500 mg/l, meets ISO 15705:2002 High range:

0 - 15000 mg/l

<sup>\*</sup> HACH® is a registered trademark of Hach Company, Loveland, Colorado. The use of the HACH® trademark does not imply any affiliation with or approval by Hach Company regarding the formulation, testing or compatibility of these products for use in HACH® brand spectrophotometers or other devices or systems.

### Thermoreactor RD 125

# For the Tube test digestion of:

COD (150 °C) TOC (120 °C) Total Chromium (100 °C) Total Nitrogen (100 °C) Total Phosphate (100 °C)



Chemical digestion of samples is required for the photometric determination of COD, TOC, total phosphate and total nitrogen.

The required temperatures and reaction time can be selected by using the membrane keypad of the reactor RD 125. The unit works at three different temperatures (100 / 120 / 150 °C) and three pre-set reaction times 30 / 60 / 120 minutes). When digestion is complete, the reactor automatically switches off and gives a corresponding LED indication with short beep alarm.

The RD 125 reactor is fitted with 24 holes for 16 mm diameter vials.

With the voltage switch on the back 230 V and 115 V are selectable.

**COD Reactor RD 125** Order code: 2 41 89 40

#### Technical data RD 125

115 V / 50-60 Hz (switchable)   Power   550 W     Dimensions   248 x 219 x 171 mm     Weight   3.9 kg     Materials, housing   PPS     Lid   PC     Block insert   PBT     Heating block   Aluminium     Holes in the   24 holes,     aluminium block   16.2 mm ± 0.2 mm     Selectable temp.   100 / 120 / 150 °C     Probe type   Pt100 A class     Temperature   ± 1 °C at the Pt100     stability     Selected time   30 / 60 / 120 / min.     and continuous operation (∞)     Heating up   from 20 °C to 150 °C     in 12 min.     Protection against overheating   at the alu block     at the alu block     at 190 °C     max. 88 dB     (piezo buzzer)     Environmental   10 - 40 °C     max. 85 %		
Weight 3.9 kg  Materials, housing ABS  Protection grid PPS Lid PC  Block insert PBT  Heating block Aluminium  Holes in the 24 holes, aluminium block 16.2 mm ± 0.2 mm  Selectable temp. 100 / 120 / 150 °C  Probe type Pt100 A class  Temperature ± 1 °C at the Pt100  stability  Selected time 30 / 60 / 120 / min. and continuous operation (∞)  Heating up from 20 °C to 150 °C in 12 min.  Protection against overheating at the alu block at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental conditions 10 − 40 °C max. 85 %	Power supply Power	115 V / 50-60 Hz (switchable)
Materials, housing ABS Protection grid PPS Lid PC Block insert PBT Heating block Aluminium Holes in the aluminium block 16.2 mm ± 0.2 mm Selectable temp. 100 / 120 / 150 °C Probe type Pt100 A class Temperature ± 1 °C at the Pt100 Selected time 30 / 60 / 120 / min. and continuous operation (∞) Heating up from 20 °C to 150 °C in 12 min.  Protection against overheating at the alu block at 190 °C Beeper max. 88 dB (piezo buzzer) Environmental conditions 10 − 40 °C max. 85 %	Dimensions	248 x 219 x 171 mm
Protection grid Lid PC Block insert Heating block Holes in the 24 holes, aluminium block 16.2 mm ± 0.2 mm  Selectable temp. 100 / 120 / 150 °C Probe type Pt100 A class  Temperature ± 1 °C at the Pt100  Selected time 30 / 60 / 120 / min. and continuous operation (∞) Heating up from 20 °C to 150 °C in 12 min.  Protection against overheating at the alu block at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental conditions 10 − 40 °C max. 85 %	Weight	3.9 kg
aluminium block  16.2 mm ± 0.2 mm  Selectable temp.  100 / 120 / 150 °C  Probe type  Pt100 A class  Temperature stability  Selected time  30 / 60 / 120 / min. and continuous operation (∞)  Heating up  from 20 °C to 150 °C in 12 min.  Protection against overheating  Beeper  max. 88 dB (piezo buzzer)  Environmental conditions  10 - 40 °C max. 85 %	Materials, housing Protection grid Lid Block insert Heating block	PPS PC PBT
Probe type Pt100 A class  Temperature ± 1 °C at the Pt100  stability  Selected time 30 / 60 / 120 / min. and continuous operation (∞)  Heating up from 20 °C to 150 °C in 12 min.  Protection against overheating at the alu block at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental conditions 10 − 40 °C max. 85 %	Holes in the aluminium block	<b>,</b>
Temperature stability  Selected time 30 / 60 / 120 / min. and continuous operation (∞)  Heating up from 20 °C to 150 °C in 12 min.  Protection against overheating at the alu block at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental conditions 10 - 40 °C max. 85 %	Selectable temp.	100 / 120 / 150 °C
stability  Selected time  30 / 60 / 120 / min. and continuous operation (∞)  Heating up  from 20 °C to 150 °C in 12 min.  Protection against overheating  at the alu block at 190 °C  max. 88 dB (piezo buzzer)  Environmental conditions  30 / 60 / 120 / min. and continuous operation (∞)  from 20 °C to 150 °C in 12 min.  at the alu block at 190 °C  max. 88 dB (piezo buzzer)	Probe type	Pt100 A class
and continuous operation (∞)  Heating up from 20 °C to 150 °C in 12 min.  Protection against overheating at the alu block at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental conditions 10 – 40 °C max. 85 %	Temperature stability	± 1 °C at the Pt100
in 12 min.  Protection against at the alu block overheating at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental conditions 10 – 40 °C max. 85 %	Selected time	
overheating at 190 °C  Beeper max. 88 dB (piezo buzzer)  Environmental 10 – 40 °C conditions max. 85 %	Heating up	
(piezo buzzer)  Environmental 10 – 40 °C conditions max. 85 %	Protection against overheating	
conditions max. 85 %	Beeper	
rei. Humaity	Environmental conditions	

**CE-Conformity** 

## Waste Water Setups

**Waste Water Setup MD 600** 21 41 00 Photometer MD 600 with standard accessory, Infrared data transmission module IRiM

**Waste Water Setup MD 610** 21 41 10 Photometer MD 600 with standard accessory **Bluetooth®** data transmission

Waste Water Setup SpectroDirect 71 21 0 Spectrophotometer SpectroDirect with standard accessory, 5 round vials Ø 24 mm

#### **Delivery Content**

- Thermoreactor RD 125
- tube stand
- membrane filter set
- instruction manual
- warranty information

#### Ranges

COD 0 - 150 mg/l and 0 - 1500 mg/l, Ammonia 1 - 50 mg/l N, Nitrate 1 - 30 mg/l N Nitrite LR 0,01 - 0,3 mg/l N Nitrogen 5 - 150 mg/l N Phosphate 0.02 - 1 mg/l P / 0.06 - 3.5 mg/l PO<sub>4</sub>

Reagents	
COD 0-150 mg/l O <sub>2</sub> (25 pc.), mercury free ** (25 pc.) (150 pc.)	2 42 07 10 2 42 07 20 2 42 07 25
CSB 0-1500 mg/l $O_2$ (25 pc.), mercury free ** (150 pc.), mercury free ** (25 pc.) (150 pc.)	2 42 07 11 2 42 07 16 2 42 07 21 2 42 07 26
CSB 0-15000 mg/l $O_2$ (25 pc.), mercury free ** (25 pc.) (150 pc.) ** without chloride removal	2 42 07 12 2 42 07 22 2 42 07 27
Ammonia VARIO HR tube test	53 56 50
Nitrate VARIO tube test	53 55 80
Nitrite LR VARIO powder pack	53 09 80
Nitrogen VARIO Total HR tube test	53 55 60
Phosphate VARIO Total HR tube test	53 52 10

Accessories	
Set of round vials with lids Height 48 mm, Ø 24 mm	19 76 29
Membrane filter set for use when preparing samples, 25 membra filters 0.45 μm, 2 syringes 20 ml	36 61 50 ane
Vial stand for 6 round vials Ø 24 mm, acrylic glass	41 89 51
Vial stand for 10 vials (Ø 16 mm or □ 13,5 mm), acrylic gla	41 89 57 ss
Automatic pipette*, 1 - 5 ml	41 90 76
Pipette tips*, 1 - 5 ml (white), 100 pc.	41 90 66
Automatic pipette**, 0.1 - 1 ml	41 90 77
Pipette tips**, 0,1 - 1 ml (white), 1000 pc.	41 90 73
* 0 - 150 mg/l and 0 - 1500mg/l ; ** 0 - 15000 mg/l	