Photometer MD 600 & MD 610



Highlights

- Highest/reproducible precision with interference filter
- Display with background lighting
- More than 120 pre-programmed methods
- Automatic selection of wavelength
- User guidance in German, English,
 French, Spanish, Italian, Portuguese (BR),
 Polish, and Indonesian.
- Buffer for 1000 data records (MD 600),
 500 data records (MD 610)

- More than 35 user-specific methods possible
- Bluetooth® interface for connection to smart phones and tablets (only with MD 610)*
- iOS® and Android™ app for data management and email delivery (only with MD 610)*
- Infrared interface (only with MD 600)
- Waterproof housing*
- Handheld format, portable

www.lovibond.com

68

^{*)} analog IP 68, 1 Stunde bei 0,1 Meter

The MD 610 and MD 600 give you mobile devices in a modern design with the analytical features of laboratory photometers.

All important water analysis parameters from A(luminium) to Z(inc) are covered by these two devices. Combined with the high precision of Lovibond® reagents, a reliable and quick analysis of water samples is guaranteed. Reagent tablets, powder reagents, liquid reagents, or cuvette tests are used depending on the method.

Six long-lasting LEDs serving as a light source in combination with interference filters guarantee the highest precision. The devices are designed without moving optical parts and thus have a maintenance-free measuring unit. Up to 1,000 data records can be stored in the MD 600 (500 data records in the MD 610).

The **AqualX®** app, available free of charge, offers the possibility of transferring measurements to smart phones or tablets via **Bluetooth®**. The data management then enables analysis and export as a CSV file or graph via email. The app is available free of charge for Android™ and iOS®.

The proven MD 600 photometer uses the classic infrared interface with which data can be transferred by means of the IRIM module to the PC or laptop.



The instrument has a factory calibration, which is related to international standards, which are not N.I.S.T traceable. The instrument may be calibrated by the user in a "user calibration mode" with N.I.S.T traceable standards.

(N.I.S.T. = National Institute of Standards and Technology)

New methods

Test methods are regularly updated to suit market requirements. You can find software updates for new methods and additional languages on our website at www.lovibond.com.

Polynomials

With the help of an external mathematical program, the corresponding polynomial is created from data pairs (concentration/absorption). A known polynomial may also be used. 25 order polynomials ($y = A + Bx + Cx^2 + Dx^3 + EX^4 + FX^5$) can be stored together with user-specific parameters such as wavelength, measuring range, unit and number of decimals.

Concentration

This function can be used to measure 2 to 14 known standards. On the basis of the concentrations/absorption pairs obtained, the photometer will calculate a linear interpolation between the measured points. Up to 10 methods can be stored for further sample measurements.

Delivery Content

- Instrument in carrying case
- 4 batteries
- 3 Round vials each 24 and 16 mm ø
- 1 adapter each for 16 mm and 13 mm vials
- Plastic stirring rod 13 cm,
 Brush 11 cm, screw driver
- Warranty information
- Certificate of Compliance
- Instruction Manual

Order codes (without reagents)

MD 600: 21 40 20 MD 610: 21 40 25

Please specify the reagents or parameters required at time of order.

You can find updated information on parameters and measuring ranges on our website at www.lovibond.com

Applications

- Waste Water
- Drinking Water
- Industrial Process Water
- Science & Research
- Governmental and Private Laboratories
- Mobile Applications



Verification Standard Kit

The verification standard kit for the MD 600 / 610 is designed to assure the user of the accuracy and the reliability of the results related to the integrated wave lengths.

The shelf life of the verification standard kit is two years from the date of production, provided that storage and use are in accordance with the instructions provided.

Verification Standard Kit

21 56 40

0

Please see pages 78 onwards for reagents (order codes)



Photometer MD 600 & MD 610



Technical Data			Accessories		
Display	Backlit graphic-display	Operation	Acid and solvent resistant, touch-sensitive keypad with audibe feedback via integrated	Item	Code
Interfaces	Infrared¹ (MD 600), Bluetooth® 4.0 (MD 610) RJ45 socket for Internet updates²			Set of 12 round vials with lid Height 48 mm, Ø 24 mm	19 76 20
		Power Supply	beeper 4 batteries (Mignon AA/LR6); Operation time: approx. 26 h continuous operation or 3500 tests	Set of 10 round vials with lid Height 90 mm, Ø 16 mm	19 76 65
Optics	LEDs, interference filters (IF) and photo sensor in transparent sample chamber Wavelength range: 430 nm IF $\Delta \lambda = 5 \text{ nm}$ 530 nm IF $\Delta \lambda = 5 \text{ nm}$ 560 nm IF $\Delta \lambda = 5 \text{ nm}$ 580 nm IF $\Delta \lambda = 5 \text{ nm}$ 610 nm IF $\Delta \lambda = 6 \text{ nm}$ 660 nm IF $\Delta \lambda = 5 \text{ nm}$ IF = interference filter			Adapter for round vials ø 16 mm	19 80 21 90
				Adapter for round vials ø 13 mm	19 80 21 92
		Auto-Off	approx. 20 minutes after last keypress with audible signal	Set of multy vials-3 with lids path length 10 mm, 10 ml volume Height 48 mm, Ø 24 mm (12 pc.)	19 76 05
		Dimensions	approx. 210 x 95 x 45 mm (unit) approx. 395 x 295 x 106 mm (case)	Vial stand for 6 round vials Ø 24 mm, acrylic glass	41 89 51
		Weight (unit)	approx. 450 g	Vial stand for 10 vials (Ø 16 mm or □ 13,5 mm), acrylic gla	41 89 57 ss
		Ambient Conditions	5–40 °C at max. 30–90 % rel. humidity (non condensing)	Sealing ring for vial ø 24 mm (12 pc.)	19 76 26
Wavelength Accuracy	± 1 nm			Battery, 1.5 V, AA-Alkali-Mangan (4 pc.)) 19 50 025
		Language		Cleaning cloth for vials	19 76 35
Photometric Accuracy*	2 % FS (T = 20 °C – 25 °C)	Language Selection	German, English, French, Spanish, Italian,Portuguese, Polish, Indonesian ; additional languages via Internet update	Plastic funnel with handle	47 10 07
				Plastic stirring rod, 13 cm length	36 41 00
Photometric Resolution	0.005 A			Plastic stirring rod, 13 cm length, (10 pc.)	36 41 20
		Memory Capacity	approx. 1000 data sets (MD 600) approx. 500 data sets (MD 610)	Plastic stirring rod, 10 cm length	36 41 09
				Plastic stirring rod, 10 cm length, (10 pc.)	36 41 30
		CE-Conformity		Cleaning brush, 10 cm	38 02 30
		1 antional auxilable. IDIA / /afrayad laterfore Nandul)		Verification Standard Kit	21 56 40
Please see pages 78 onwards for reagents (order codes)		 optional available: IRiM (Infrared Interface Modul) optional available: connection cable with integrated electronics (RS 232 / RJ-45 plug) tested with standard solutions 		Cable for update for connection to a PC	21 40 30
				Data transmission modul IRiM	21 40 50

The **Bluetooth®** word mark is a registered trademark owned by Bluetooth SIG, Inc. and any use by Lovibond® Tintometer GmbH is under license. IOS® is a registered trademark of Cisco, Inc. and licensed to Apple, Inc. Android™ is a trademark of Google, Inc.

70



Infrared data transmission modul IRiM



The IRiM (infrared interface modul) uses modern infrared technology to transmit measurement data from the MD 600 photometer to one of 3 optional interfaces. These interfaces can be used to connect to a PC, a USB printer¹⁾ or alternatively a serial printer²⁾. The interface which is selected is displayed by an LED function indicator. The user can switch between the interfaces using the "Select" button.

The unit is supplied complete with data logging software providing easy and rapid transfer of data to the PC. As an option, the data can be saved as an Excel sheet or a .txt file.

Measurement data can quickly be printed out, using a specified¹⁾ USB or alternatively a printer with a serial plug-in connected to the IRiM.

Applicable for the following operating systems: Windows XP, Windows Vista and Windows 7/10.

1) USB printer: HP Deskjet 6940; 2) each ASCII printer

Delivery content

The IRiM is delivered ready for use, with the following accessories :

USB cable, 4 batteries, screwdriver, CD-ROM, operating instructions and guarantee certificate

Order code: 21 40 50

Photometer and Fluorometer MD 640



Highlights

- Inbuilt PTSA & Fluorescein measurement - no adapter required
- High quality results due to interference filters and long-life LEDs
- Automatic wavelength selection
- Bluetooth® data transmittance to Lovibond® AquaLX® App
- Covers more than 120 important methods for water analysis such as aluminium, chlorine, COD, bromine, chlorine dioxide, copper, iron, molybdate and phosphate
- Advanced data management via AquaLX® App
- Portable and easy handling
- One time zero
- Data storage for 500 data sets
- Robust, water proof design
- Backlit display

Introduction The Lovibond® Photometer MD 640 is an

enhanced version of the MD 610 photometer, offering additional fluorescence capability for the determination of PTSA and fluorescein in water systems.

PTSA (1,3,6,8 pyrenetetrasulfonic acid, sodium salt) and fluorescein are fluorescent materials that are increasingly being added to speciality water treatment products to enable real time product dose analysis. Both materials are detectable at ppb levels, are non-toxic and chemically stable, all of which make them ideal tracer additives throughout complex water systems. Accurately measuring product dose levels helps the water treatment specialist to control water chemistry; prevent corrosion, scale and biological fouling; increase system efficiency and, ultimately, save energy and costs.

Delivery Content

- Instrument in carrying case
- 4 hatteries
- 3 round vials each 24 and 16 mm ø (black lid)
- 1 adapter each for 16 mm and 13 mm vials
- Plastic stirring rod 13 cm, Brush 11 cm, syringe 5 ml, screw driver
- Warranty information
- Certificate of Compliance
- Instruction Manual

Order codes (without reagents) MD 640: 21 41 40

Please specify the reagents or parameters required at time of order.

You can find updated information on parameters and measuring ranges at www.lovibond.com

Applications

- Industrial Process Water & Waste Water
- **Drinking Water**
- Science & Research
- Governmental and Private Laboratories
- Mobile Applications

Technical Data

Bluetooth® 4.0 **Interfaces** RJ45 socket for Internet updates¹

Optics

Display

photo sensor in transparent sample chamber Wavelength range: 430 nm IF $\Delta \lambda = 5$ nm 530 nm IF $\Delta \lambda = 5$ nm 560 nm IF $\Delta \lambda = 5$ nm 580 nm IF $\Delta \lambda = 5$ nm 610 nm IF $\Delta \lambda = 6$ nm

660 nm IF $\Delta \lambda = 5$ nm IF = interference filter

LEDs, interference filters (IF) and

Backlit graphic-display

UV excitation 375 nm

PTSA 10 - 1000 ppb Measurement Ranges Fluorescein 10 - 400 ppb

Calibration Monthly (user) Check (using calibration sets)

Calibration Factory set & user adjustable (using calibration Standard Set)

Wavelength Accuracy

Photometric 2 % FS (T = 20 °C - 25 °C)Accuracy*

± 1 nm

Photometric

0.005 A Resolution

Operation

Acid and solvent resistant, touch-sensitive keypad with audible feedback via integrated

Power Supply

4 batteries (Mignon AA/LR6); Operation time: approx. 26 h continuous operation or 3500 tests

Auto-Off

Approx. 20 minutes after last keypress with audible signal Approx. 210 x 95 x 45 mm (unit)

Dimensions

approx. 395 x 295 x 106 mm (case) Approx. 450 g

Weight (unit) Ambient

5-40 °C at max. 30-90 % Conditions rel. humidity (non condensing)

Language Selection

German, English, French, Spanish, Italian, Portuguese, Polish, Indonesian; additional languages via Internet update

Memory Capacity

Approx. 500 data sets

CE-Conformity

Code ltem Set of 12 round vials with lid 19 76 20 Height 48 mm, Ø 24 mm Set of 12 round vials with black lid 19 76 57 for PTSA / Fluorescein Height 48 mm, Ø 24 mm Set of 10 round vials with lid 19 76 65 Height 90 mm, Ø 16 mm Adapter for round vials ø 16 mm 19 80 21 90 Adapter for round vials ø 13 mm 19 80 21 92 Set of **multi vials-3** with lids 19 76 05 path length 10 mm, 10 ml volume Height 48 mm, Ø 24 mm (12 pc.) Vial stand for 6 round vials 41 89 51 Ø 24 mm, acrylic glass Vial stand for 10 vials 41 89 57 (Ø 16 mm or □ 13,5 mm), acrylic glass Sealing ring for vial ø 24 mm (12 pc.) 19 76 26 Battery, 1.5 V, AA-Alkali-Mangan (4 pc.) 19 50 025 Cleaning cloth for vials 19 76 35 Plastic funnel with handle 47 10 07 Plastic stirring rod, 13 cm length 36 41 00 Plastic stirring rod, 13 cm length, (10 pc.) 36 41 20 Plastic stirring rod, 10 cm length 36 41 09 Plastic stirring rod, 10 cm length, (10 pc.) 36 41 30 Cleaning brush, 10 cm 38 02 30 Verification Standard Kit 21 56 40 Cable for update 21 40 30 for connection to a PC PTSA standard addition solution, 46 12 10 1000 ppb, 50ml PTSA calibration set 46 12 45 (0, 200, 1000 ppb) Fluorescein standard addition 46 12 30 solution, 400 ppb, 50ml Fluorescein calibration set 46 12 40 (0, 75, 400 ppb)

Accessories

¹ optional available: connection cable with integrated electronics (RS 232 / RJ-45 plug)

^{*} tested with standard solutions