

# Dip & Read Tests

For the Dairy Industry



## Rapid Tests for the Dairy Industry

- Fast results
- Easy to use
- Precise readings



**MACHERY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)








Since 1911

# MN Dip & Read Tests for the Dairy Industry



The dairy industry uses test strips as easy-to-use tools mainly in quality control. The Dip & Read tests in this special selection provide valuable information for dairy farms and the milk-product industry. Test strips for peroxide or peracetic acid are used to ensure that sanitizers have been fully rinsed from packages prior to filling. Heating processes of milk are easily checked using test strips for the enzymes phosphatase and peroxidase.

## Product information


Test	Product	Application	Gradation	REF
	QUANTOFIX® Peroxide 25	<ul style="list-style-type: none"> <li>• Detects peroxide traces down to 0.5 mg/l</li> <li>• Strips can be stored at room temperature</li> <li>• Protected against false positive results</li> </ul>	0 · 0.5 · 2 · 5 · 10 · 25 mg/l H <sub>2</sub> O <sub>2</sub>	<b>913 19</b>
	QUANTOFIX® Peracetic acid 50	<ul style="list-style-type: none"> <li>• Specific for peracetic acid</li> <li>• No peroxide interference during reaction time</li> <li>• Strips can be stored at room temperature</li> </ul>	0 · 5 · 10 · 20 · 30 · 50 mg/l peracetic acid	<b>913 40</b>
	Phosphatesmo MI	<ul style="list-style-type: none"> <li>• Detection of alkaline phosphatase in milk</li> <li>• Easy quality control for pasteurization</li> <li>• Easy Dip &amp; Read procedure</li> </ul>	Detection limit: 300 U/l alkaline phosphatase in UHT milk	<b>906 12</b>
	Peroxtesmo MI	<ul style="list-style-type: none"> <li>• Detection of lactoperoxidase in milk</li> <li>• Easy quality control for high heated milk</li> <li>• Odorless and non-toxic for maximum usability</li> </ul>	Detection limit: 3% of raw milk in UHT milk	<b>906 27</b>
	pH-Fix 6.0-7.7	<ul style="list-style-type: none"> <li>• Easy Dip &amp; Read test on milk pH</li> <li>• 12 other pH ranges also available</li> <li>• Fixed indicator does not contaminate sample</li> </ul>	pH 6.0 · 6.4 · 6.7 · 7.0 · 7.3 · 7.7	<b>921 50</b>

### Others

WATOR	Water distribution in butter for easy quality assessment	<b>906 10</b>
Udder test paper	Easy indicator for mastitis	<b>907 48</b>
QUANTOFIX® Peroxide 1000	Detects peroxide levels up to 1000 mg/l	<b>913 33</b>
QUANTOFIX® Calcium	Easy detection of 10-100 mg/l calcium in milk	<b>913 24</b>
QUANTOFIX® Chlorine Sensitive	Detection of total chlorine traces down to 0.1 mg/l	<b>913 39</b>

## Contact

Your local distributor



# MACHEREY-NAGEL

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany  
 Germany and international: Switzerland: France: USA  
 Tel.: +49 (0) 24 21 96 90 Tel.: +41 (0) 62 388 55 00 Tel.: +33 (0) 3 88 68 22 68 Tel.: +1 484 821 0984  
 Fax: +49 (0) 24 21 96 91 99 Fax: +41 (0) 62 388 55 05 Fax: +33 (0) 3 88 51 76 88 Fax: +1 484 821 1272  
 e-mail: sales-de@mn-net.com e-mail: sales-ch@mn-net.com e-mail: sales-fr@mn-net.com e-mail: sales-us@mn-net.com

