

# Dip & Read Tests

For the Beverage Industry



## Rapid Tests for the Beverage Industry

- Fast & easy
- Reliable results
- Efficient quality control



**MACHEREY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)









Since 1911

# MN Dip & Read Tests for the Beverage Industry

For beverage manufactures and breweries reliable disinfection and sterility of surfaces, machinery, bottles and containers is vital for product quality. For disinfection, often chemical sanitizer solutions are used, which have to be carefully rinsed out so they don't contaminate the product. Test strips provide easy means to check the disinfection and rinsing process. Additionally, dip and read tests are often used for Quality control. MN test strips are ideal for many applications in the beverage industry and already provide companies around the globe with quick and reliable results.



## Product information

Test	Product	Application	Gradation	REF
	QUANTOFIX® Peroxide 25	<ul style="list-style-type: none"> <li>• Detects peroxide residues down to 0.5 mg/l</li> <li>• Strips can be stored at room temperature</li> <li>• Protection against false positive results</li> </ul>	0 · 0.5 · 2 · 5 · 10 · 25 mg/l H <sub>2</sub> O <sub>2</sub>	<b>913 19</b>
	QUANTOFIX® Peracetic acid 50	<ul style="list-style-type: none"> <li>• Specific to peracetic acid</li> <li>• No peroxide interference during reaction time</li> <li>• Strips can be stored at room temperature</li> </ul>	0 · 5 · 10 · 20 · 30 · 50 mg/l peracetic acid	<b>913 40</b>
	QUANTOFIX® QUAT	<ul style="list-style-type: none"> <li>• Easy detection of potency levels</li> <li>• Conversion table for many quaternary ammonium compounds (QUATs) provided</li> <li>• Ensures optimal disinfection</li> </ul>	0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/l benzalkonium chloride	<b>913 37</b>
	QUANTOFIX® Peroxide 1000	<ul style="list-style-type: none"> <li>• Detection of peroxide potency levels up to 1000 mg/l</li> <li>• Second test pad indicates concentration &gt;2000 mg/l</li> <li>• Ensures optimal disinfection</li> </ul>	0 · 50 · 150 · 300 · 500 · 800 · 1000 mg/l H <sub>2</sub> O <sub>2</sub>	<b>913 33</b>
	QUANTOFIX® Peracetic acid 2000	<ul style="list-style-type: none"> <li>• Specific to peracetic acid</li> <li>• No peroxide interference during reaction time</li> <li>• Detects potency levels up to 2000 mg/l</li> </ul>	0 · 500 · 1000 · 1500 · 2000 mg/l peracetic acid	<b>913 42</b>
	QUANTOFIX® Ascorbic acid	<ul style="list-style-type: none"> <li>• Easy detection of vitamin C in juice</li> <li>• Easy Dip &amp; Read procedure provides results in 30 seconds</li> </ul>	0 · 50 · 100 · 200 · 300 · 500 · 1000 · 2000 mg/l vitamin C	<b>913 14</b>
<b>Others</b>				
	QUANTOFIX® Cyanide	Cyanide control during production of fruit brandy		<b>913 18</b>
	QUANTOFIX® Sulphite	Control of sulphite in wine and foodstuff treated with sulphur compounds for longer shelf life		<b>913 06</b>
	QUANTOFIX® Nitrate/Nitrite	Detection of nitrate (10-500 mg/l) and nitrite (1-80 mg/l)		<b>913 13</b>
	QUANTOFIX® Tin	Detection of tin (10-500 mg/l) in canned beverages		<b>913 09</b>

## Contact

Your local distributor:



# MACHEREY-NAGEL

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany  
 Germany and international: Tel.: +49 (0) 24 21 96 90 Fax: +49 (0) 24 21 96 91 99 e-mail: sales-de@mn-net.com  
 Switzerland: Tel.: +41 (0) 62 388 55 00 Fax: +41 (0) 62 388 55 05 e-mail: sales-ch@mn-net.com  
 France: Tel.: +33 (0) 3 88 68 22 68 Fax: +33 (0) 3 88 51 76 88 e-mail: sales-fr@mn-net.com  
 USA: Tel.: +1 484 821 0984 Fax: +1 484 821 1272 e-mail: sales-us@mn-net.com

