

Dip & Read Tests

For Healthcare



Rapid Tests for Healthcare

- Quick & easy
- Reliable results
- Higher patient safety



MACHEREY-NAGEL

www.mn-net.com








Since 1911

MN Dip & Read Tests for Healthcare

Health care professionals use different types of disinfectants and sanitizers for cleaning of medical devices. After disinfection all devices have to be rinsed and controlled for sanitizer residues prior to the next use. Dip & Read tests allow the quick and easy detection of sanitizer residues. Thus, MN test strips increase patient safety, while making healthcare professionals' work easier and more efficient.



Product information

Test	Product	Application	Gradation	REF
	QUANTOFIX® Peroxide 100 CE	<ul style="list-style-type: none"> • Easy detection of peroxide residues • CE-marked (medical device directive 93/42/EWG) • Strips can be stored at room temperature 	0 · 1 · 3 · 10 · 30 · 100 mg/l H ₂ O ₂	913 12
	QUANTOFIX® Peracetic acid 50	<ul style="list-style-type: none"> • Specific for peracetic acid • No peroxide interference during reaction time • Ideal for peracetic acid based sanitizers 	0 · 5 · 10 · 20 · 30 · 50 mg/l peracetic acid	913 40
	QUANTOFIX® QUAT	<ul style="list-style-type: none"> • Easy detection of potency levels of sanitizers • Conversion table for many quaternary ammonium compounds (QUATs) provided • Ensures optimal disinfection 	0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/l benzalkonium chloride	913 37
	pH-Fix 3.6-6.1 CE	<ul style="list-style-type: none"> • Easy control of rinsing after disinfection with hydroxyacetic acid • CE-marked (medical device directive 93/42/EWG) • Fixed indicator does not contaminate sample 	pH 3.6 · 4.1 · 4.4 · 4.7 · 5.0 · 5.3 · 5.6 · 6.1	921 50
	QUANTOFIX® Chlorine Sensitive CE	<ul style="list-style-type: none"> • Easy detection of residues of chlorine in feed water and rinse water • CE-marked (medical device directive 93/42/EWG) 	0 · 0.1 · 0.5 · 1 · 3 · 10 mg/l Cl ₂	913 39
Others	Medi-Test Urinalysis	A comprehensive range of urine dipsticks and URYXXON® instruments for use in physician offices as well as hospitals		930 ...

Contact

Your local distributor



MACHERY-NAGEL

MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany
 Germany and international: Tel.: +49 (0) 24 21 96 90
 Switzerland: Tel.: +41 (0) 62 388 55 00
 France: Tel.: +33 (0) 3 88 68 22 68
 USA: Tel.: +1 484 821 0984
 Fax: +49 (0) 24 21 96 91 99
 Fax: +41 (0) 62 388 55 05
 Fax: +33 (0) 3 88 51 76 88
 Fax: +1 484 821 1272
 e-mail: sales-de@mn-net.com
 e-mail: sales-ch@mn-net.com
 e-mail: sales-fr@mn-net.com
 e-mail: sales-us@mn-net.com



Image Credits © Monkey Business, vale, mario beauregard, - Fotolia.com
 KATEN100096/Healthcare en/12/10/2009 PD Printed in Germany