


Dip & Read Tests

For the Food Industry



Rapid Tests for the Food Industry

- Fast & easy
- Reliable results
- Efficient quality control



MACHEREY-NAGEL

www.mn-net.com



Since 1911

MN Dip & Read Tests for the Food Industry

High quality and consumer safety is a fundamental aspect in food production. MN offers a wide range of test strips that can be used to verify constant high quality of products and safety of equipment and production facilities. MN test strip solutions are fast, easy to use and reliable, making products and operations safe and efficient at the same time.



Product information

Test	Product	Application	Gradation	REF
	QUANTOFIX® Peroxide 25	- Reliable detection of peroxide residues down to 0.5 mg/l - Can be stored at room temperature - Accurate results in 60 seconds	0 · 0.5 · 2 · 5 · 10 · 25 mg/l H ₂ O ₂	913 19
	QUANTOFIX® Sulphite	- Sulphite in food can cause severe allergic reactions - Tests for sulphite in food treated with sulphur based preservatives - Reliable results in 20 seconds	0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/l SO ₃ ²⁻	913 06
	QUANTOFIX® Peracetic Acid 2000	- Specific to peracetic acid - No peroxide interference during reaction time - Detects disinfectant potency levels up to 2000mg/l	0 · 500 · 1000 · 1500 · 2000 mg/l peracetic acid	913 42
	pH-Fix 0-14	- Wide range of applications - Bright colour scale - Reliable results in seconds	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	921 10
	QUANTOFIX® Ascorbic Acid (Vitamin C)	- Vitamin C is important and desirable in food - Also used as stabilizer and reducing agent - Reliable results in 30 seconds	0 · 50 · 100 · 200 · 300 · 500 · 1000 · 2000 mg/l vitamin C	913 14
	QUANTOFIX® Nitrate/Nitrite	- Prevalent in many foods - High doses have adverse health effects - Reliable results in 60 seconds	0 · 10 · 25 · 50 · 100 · 250 · 500 mg/l NO ₃ ⁻ 0 · 1 · 5 · 10 · 20 · 40 · 80 mg/l NO ₂ ⁻	913 13
Others				
	Peroxtesmo KO	Quick and easy detection of peroxidase in food		906 06
	PEHANON	Easy pH determination in coloured liquids		904 ...
	QUANTOFIX® Cyanide	Reliable cyanide detection, especially in fruit brandy production		913 18
	INDIPRO	Detection of protein residue on equipment and working surfaces		907 65

Contact

Your local distributor:



MACHERY-NAGEL

MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany
 Germany and international: Tel.: +49 (0) 24 21 96 90 Fax: +49 (0) 24 21 96 91 99 e-mail: sales-de@mn-net.com
 Switzerland: Tel.: +41 (0) 62 388 55 00 Fax: +41 (0) 62 388 55 05 e-mail: sales-ch@mn-net.com
 France: Tel.: +33 (0) 3 88 68 22 68 Fax: +33 (0) 3 88 51 76 88 e-mail: sales-fr@mn-net.com
 USA: Tel.: +1 484 821 0984 Fax: +1 484 821 1272 e-mail: sales-us@mn-net.com

